

Pill Dispenser Monitor

The first button – turn on/off “wrong pill” mode.

The second button – resets and clears status and pill flags.

Temperature gauges between (15-40).


RTC – Real time clock settings.

LED Status – Indicates 4 different device status:

1. Blue – neutral status, temperature ok, no pill taken.
2. Green – temperature ok, correct pill has been taken.
3. Yellow – temperature ok, wrong pill has been taken.
4. Red – temperature is not ok, dispenser overheating.

2xWebhook:

1. To send temperature information to the Integromat.
2. To send pill status information to the Integromat.



Gets the values -> logs them in spreadsheets based on the case -> notify the telegram group with the correct message.